

Risedronate Sodium

Catalog No: #S1041


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Risedronate Sodium |
| Brief Description | Inhibitors |
| Purification | 99.83% |
| Target Name | HBV inhibitor; Others inhibitor |
| Calculated MW | 305.09 |
| Formulation | C7H10NNaO7P2 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Metabolism

system SMILES: O.O.O.O.O.[Na+].OC(CC1=CN=CC=C1)(P(O)(O)=O)P(O)([O-])=O.OC(CC1=CN=CC=C1)(P(O)(O)=O)P(O)([O-])=O

Pathways: Microbiology & virology; Others Receptor: HAP; FPS Boling pt: Melting pt: 692.3C Solubility: Water: 2 mg/mL (6.55 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Bergstrom JD, et al. Arch Biochem Biophys. 2000 Jan 1;373(1):231-41.

Note: This product is for in vitro research use only