

## Salifungin


Catalog No: #S1043

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Salifungin
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	326.57
Formulation	C13H9BrClNO2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: OC1=C(C=C(Br)C=C1)C(=O)NC1=CC=C(Cl)C=C1  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: Solubility: DMSO: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Kulagin VI, et al. Vestn Dermatol Venerol. 1972 Aug;46(8):75-6.

Note: This product is for in vitro research use only