

Clopidol

Catalog No: #S1044


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Clopidol
Brief Description	Inhibitors
Purification	99.00%
Target Name	Others
Calculated MW	192.04
Formulation	C7H7Cl2NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC1=C(Cl)C(=O)C(Cl)=C(C)N1 Pathways: Others Receptor: Others Boiling pt: Melting pt: >320°C Solubility: DMSO: 1 mg/mL (5.2 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Fang B, et al. Anal Sci, 2009, 25(10), 1203-1206.

Note: This product is for in vitro research use only