

Proglumide


Catalog No: #S1052

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Proglumide
Brief Description	Inhibitors
Purification	99.44%
Target Name	cholecystokinin antagonist
Calculated MW	334.42
Formulation	C18H26N2O4
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
 SMILES: CCCN(CCC)C(=O)C(CCC(O)=O)NC(=O)C1=CC=CC=C1
 Pathways: GPCR/G Protein Receptor: CCKB
 Boiling pt: Melting pt: 148-150°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Bunney BS, et al. Ann N Y Acad Sci. 1985;448:345-51.

Note: This product is for in vitro research use only