

Mebeverine hydrochloride

Catalog No: #S1053


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Mebeverine hydrochloride |
| Brief Description | Inhibitors |
| Purification | 99.00% |
| Target Name | Others |
| Calculated MW | 466.01 |
| Formulation | C25H36ClNO5 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Other
 SMILES: Cl.CCN(CCCOC(=O)C1=CC(OC)=C(OC)C=C1)C(C)CC1=CC=C(OC)C=C1
 Pathways: Others
 Receptor: Others
 Boiling pt: Melting pt: 105-107°C
 Solubility: DMSO: 10 mM
 Appearance: White Powder
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Rahman MZ, et al. Mymensingh Med J. 2014 Jan;23(1):105-13.

Note: This product is for in vitro research use only