

Azilsartan

Catalog No: #S1057


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Azilsartan
Brief Description	Inhibitors
Purification	98.00%
Target Name	RAAS antagonist
Calculated MW	456.46
Formulation	C25H20N4O5
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: CCOC1=NC2=CC=CC(C(O)=O)=C2N1CC1=CC=C(C=C1)C1=CC=CC=C1C1=NC(=O)ON1 Pathways: Endocrinology/Hormones Receptor: AT1 receptor Boiling pt: Melting pt: Solubility: DMSO: 91 mg/mL (199.36 mM) Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ojima M, et al. J Pharmacol Exp Ther, 2011, 336(3), 801-808.

Note: This product is for in vitro research use only