

Retinol palmitate

Catalog No: #S1059


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Retinol palmitate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	524.86
Formulation	C36H60O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CCCCCCCCCCCCCCCC(=O)OC/C=C(\C)/C=C/C=C(\C)/C=C/C1=C(CCCC1(C)C)C Pathways: Others Receptor: Others Bolin
 g pt: 607.5°C Melting pt: 416.78°C Solubility: Soluble in hydrocarbons Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Athanassiades TJ, et al. J Natl Y Inst, 1981, 67(5), 1153-1156.

Note: This product is for in vitro research use only