

Toltrazuril

Catalog No: #S1061


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Toltrazuril
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	425.38
Formulation	C18H14F3N3O4S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CN1C(=O)NC(=O)N(C1=O)C1=CC(C)=C(OC2=CC=C(SC(F)(F)F)C=C2)C=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: 194-196°C Solubility: DMSO: 85 mg/mL (199.82 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Mehlhorn H, et al. Parasitol Res, 1988, 75(1), 64-66.

Note: This product is for in vitro research use only