

Ticarcillin disodium

Catalog No: #S1063


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Ticarcillin disodium
Brief Description	Inhibitors
Purification	90.00%
Target Name	Others
Calculated MW	428.39
Formulation	C ₁₅ H ₁₄ N ₂ Na ₂ O ₆ S ₂
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1(C(N2C(S1)C(C2=O)NC(=O)C(C3=CSC=C3)C(=O)[O-])C(=O)[O-])C.[Na+].[Na+] Pathways: Others Receptor: Others Bolin
 g pt: 768.3C Melting pt: Solubility: DMSO: 18mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Webb D, et al. Arch Intern Med. 1978 Nov;138(11):1618-20.

Note: This product is for in vitro research use only