

DGAT2L3 Polyclonal Antibody

Catalog No: #40840

Package Size: #40840-1 50ul #40840-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DGAT2L3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	DGAT2L3 Polyclonal Antibody detects endogenous levels of DGAT2L3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human DGAT2L3.
Target Name	DGAT2L3
Other Names	AWAT1; DGA2; DGAT2L3; Acyl-CoA wax alcohol acyltransferase 1; Diacylglycerol O-acyltransferase 2-like protein 3; Diacylglycerol acyltransferase 2; Long-chain-alcohol O-fatty-acyltransferase 1
Accession No.	Swiss-Prot: Q58HT5NCBI Gene ID: 158833
Uniprot	Q58HT5
GeneID	158833;
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/5000.


Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using DGAT2L3 Polyclonal Antibody

Western Blot analysis of 293T NIH-3T3 COLO cells using DGAT2L3 Polyclonal Antibody



Note: This product is for in vitro research use only