

Baclofen


Catalog No: #S1065

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-------------------|---|
| Product Name | Baclofen |
| Brief Description | Inhibitors |
| Purification | 98.00% |
| Target Name | GABA Receptor agonist |
| Calculated MW | 213.66 |
| Formulation | C10H12ClNO2 |
| Storage | 3 years -20°C powder; 2 years -80°C in solvent; |

Images



Product Description

Research Area: Nervous system SMILES: NCC(CC(O)=O)C1=CC=C(Cl)C=C1 Pathways: Membrane transporter/Ion channel;

Neuroscience Receptor: GABAR Boling pt: Melting pt: 208-210°C Solubility: Soluble in Water Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Garcia-Gil L, et al. Neurotoxicol Teratol. 1999 May-Jun;21(3):277-83.

Note: This product is for in vitro research use only