

## Ketanserin

Catalog No: #S1066


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Ketanserin
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor antagonist
Calculated MW	395.43
Formulation	C22H22FN3O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: FC1=CC=C(C=C1)C(=O)C1CCN(CCN2C(=O)NC3=CC=CC=C3C2=O)CC1  
 Pathways: GPCR/G Protein;  
 Neuroscience Receptor: 5-HT2A; 5-HT2C (Human); 5-HT2C (Rat)  
 Boling pt: Melting pt: 227-235°C  
 Solubility: DMSO: 2 mg/mL (5.05 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Van Nueten JM, et al. J Pharmacol Exp Ther, 1981, 218(1), 217-230.

Note: This product is for in vitro research use only