

Anthralin

Catalog No: #S1068


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Anthralin
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	226.23
Formulation	C ₁₄ H ₁₀ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: OC1=CC=CC2=C1C(=O)C1=C(C2)C=CC=C1O Pathways: Others Receptor: Others Boiling pt: 464.1°C Melting pt: 178-181°C Solubility: Chloroform: 41mg/mL Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Durdu M, et al. J Am Acad Dermatol. 2015 Apr;72(4):640-50.

Note: This product is for in vitro research use only