

Cefmetazole sodium

Catalog No: #S1070


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cefmetazole sodium
Brief Description	Inhibitors
Purification	92.05%
Target Name	Others
Calculated MW	493.52
Formulation	C ₁₅ H ₁₆ N ₇ NaO ₅ S ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CN1C(=NN=N1)SCC2=C(N3[C@@H]([C@@](C3=O)(NC(=O)CSCC#N)OC)SC2)C(=O)[O-].[Na+] Pathways: Others Receptor: Others Boiling pt: Melting pt: 2-8°C Solubility: DMSO: 39mg/mL Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Inagaki K, et al. J Parenter Sci Technol. 1993 Jan-Feb;47(1):35-9.

Note: This product is for in vitro research use only