

DL-Adrenaline Hydrochloride

Catalog No: #S1071L


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DL-Adrenaline Hydrochloride
Brief Description	Inhibitors
Purification	99.00%
Target Name	Adrenergic Receptor
Calculated MW	219.67
Formulation	C ₉ H ₁₃ NO ₃ HCl
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine system SMILES:Cl.CNCC(O)C1=CC(O)=C(O)C=C1 Pathways:Neuroscience; GPCR/G Protein Receptor:Adrenergic Receptor Boiling pt:413.1°C Melting pt:157°C Solubility:Soluble in Water Appearance:White Powder Remark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Leineweber K, et al. Naunyn Schmiedebergs Arch Pharmacol. 2007 Jul;375(5):303-9.

Note: This product is for in vitro research use only