


## Description

Product Name	Loxeen
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR antagonist
Calculated MW	391.52
Formulation	C20H25NOCH4O3S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
 SMILES: CS(=O)(=O)OC(CCN1CCCC1)(C1=CC=CC=C1)C1=CC=CC=C1  
 Pathways: Neuroscience  
 Receptor: mAChR  
 Boiling pt: 585.1°C  
 Melting pt: 152.5-155°C  
 Solubility: Soluble in DMSO  
 Appearance: White Crystal  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Waelbroeck M, et al. Br J Pharmacol. 1993 Jun;109(2):360-70.

Note: This product is for in vitro research use only