

Vitamin K1

Catalog No: #S1075


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Vitamin K1
Brief Description	Inhibitors
Purification	99.53%
Target Name	Others
Calculated MW	450.71
Formulation	C31H46O2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Metabolism systemSMILES:CC1=C(C(=O)C2=CC=CC=C2C1=O)C/C=C(\C)/CCC[C@H](C)CCC[C@H](C)CCCC(C)C

Pathways:OthersReceptor:GGCX; OsteocalcinBoiling pt:140CMelting pt:-20CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Morris DP, et al. J Biol Chem. 1993 Apr 25;268(12):8735-42.

Note: This product is for in vitro research use only