

DHHC-9 Polyclonal Antibody

Catalog No: #40841

Package Size: #40841-1 50ul #40841-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	DHHC-9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu Ms Rt
Specificity	DHHC-9 Polyclonal Antibody detects endogenous levels of DHHC-9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human DHHC-9.
Target Name	DHHC-9
Other Names	ZDHHC9; CXorf11; ZDHHC10; ZNF379; ZNF380; CGI-89; Palmitoyltransferase ZDHHC9; Zinc finger DHHC domain-containing protein 9; DHHC-9; DHHC9; Zinc finger protein 379; Zinc finger protein 380
Accession No.	Swiss-Prot: Q9Y397NCBI Gene ID: 51114
Uniprot	Q9Y397
GeneID	51114;
SDS-PAGE MW	41kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of 293T cells using DHHC-9 Polyclonal Antibody

Note: This product is for in vitro research use only