

## Nafronyl oxalate

Catalog No: #S1081


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Nafronyl oxalate
Brief Description	Inhibitors
Purification	98.00%
Target Name	5-HT Receptor antagonist
Calculated MW	473.56
Formulation	C26H35NO7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Nervous system  
 SMILES: OC(=O)C(O)=O.CCN(CC)CCOC(=O)C(CC1CCCO1)CC1=CC=CC2=CC=CC=C12  
 Pathways: GPCR/G Protein; Neuroscience  
 Receptor: 5-HT  
 Boiling pt: Melting pt: Solubility: Water: 10 mM  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Kim MJ, et al. Clin Exp Nephrol. 2013 Jun;17(3):437-8.

Note: This product is for in vitro research use only