

Benzyl isothiocyanate

Catalog No: #S1082


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Benzyl isothiocyanate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	149.22
Formulation	C8H7NS
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:OtherSMILES:S=C=NCC1=CC=CC=C1Pathways:OthersReceptor:OthersBoiling pt:242-243CMelting pt:41CSolubility:DMSO: 10 mMAppearance:LiquidRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sehrawat A, Singh SV. Mol Carcinog. 2015 Apr 7.

Note: This product is for in vitro research use only