

DHRS4 Polyclonal Antibody

Catalog No: #40842

Package Size: #40842-1 50ul #40842-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	DHRS4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC
Species Reactivity	Hu
Specificity	DHRS4 Polyclonal Antibody detects endogenous levels of DHRS4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human DHRS4.
Target Name	DHRS4
Other Names	DHRS4; Dehydrogenase/reductase SDR family member 4; NADPH-dependent carbonyl reductase/NADP-retinol dehydrogenase; CR; PHCR; NADPH-dependent retinol dehydrogenase/reductase; NRDR; humNRDR; Peroxisomal short-chain alcohol dehydrogenase; PSCD
Accession No.	Swiss-Prot: Q9BTZ2NCBI Gene ID: 10901
Uniprot	Q9BTZ2
GeneID	10901;
SDS-PAGE MW	32-34kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

Images



Western Blot analysis of HepG2 cells using DHRS4 Polyclonal Antibody

Note: This product is for in vitro research use only