

Tolnaftate

Catalog No: #S1094


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Tolnaftate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Hedgehog/Smoothened inhibitor
Calculated MW	307.42
Formulation	C19H17NOS
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:CN(C(=S)OC1=CC2=CC=CC=C2C=C1)C1=CC=CC(C)=C1Pathways:GPCR/G Protein; Stem CellsReceptor:SmoothenedBoiling pt:453.4CMelting pt:110.5-111.5CSolubility:DMSO: 62 mg/mL (201.68 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Barrett-Bee K, Dixon G. Acta Biochim Pol. 1995;42(4):465-79.

Note: This product is for in vitro research use only