

Denatonium benzoate

Catalog No: #S1098


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Denatonium benzoate
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	446.58
Formulation	C28H34N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research

Area: Other SMILES: [O-]C(=O)C1=CC=CC=C1.CC[N+](CC)(CC(=O)NC1=C(C)C=CC=C1C)CC1=CC=CC=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: 164-168°C Solubility: DMSO: 10 mM Appearance: White Crystal Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Meyerhof W, et al. Chem Senses. 2010 Feb;35(2):157-70.

Note: This product is for in vitro research use only