

Formestane

Catalog No: #S1099


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Formestane
Brief Description	Inhibitors
Purification	98.00%
Target Name	Aromatase inhibitor
Calculated MW	302.41
Formulation	C19H26O3
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area:CancerSMILES:C[C@]12CCC(=O)C(=C1CC[C@@H]3[C@@H]2CC[C@]4([C@H]3CCC4=O)C)OPathways:Endocrinology/HormonesReceptor:AromataseBoiling pt:475.4CMelting pt:199-202CSolubility:DMSO: 61 mg/mL (201.71 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Wiseman LR, Goa KL. Drugs Aging. 1996 Oct;9(4):292-306.

Note: This product is for in vitro research use only