

DNA pol  $\alpha$  Polyclonal Antibody

Catalog No: #40844

Package Size: #40844-1 50ul #40844-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	DNA pol $\alpha$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	DNA pol $\alpha$ Polyclonal Antibody detects endogenous levels of DNA pol $\alpha$ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human DNA pol $\alpha$ .
Target Name	DNA pol $\alpha$
Other Names	POLA1; POLA; DNA polymerase alpha catalytic subunit; DNA polymerase alpha catalytic subunit p180
Accession No.	Swiss-Prot: P09884NCBI Gene ID: 5422
Uniprot	P09884
GeneID	5422;
SDS-PAGE MW	165kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images

Western Blot analysis of Jurkat cells using DNA pol  $\alpha$  Polyclonal Antibody

---

Note: This product is for in vitro research use only