

## Sarafloxacin hydrochloride

Catalog No: #S1107


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sarafloxacin hydrochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	421.83
Formulation	C <sub>20</sub> H <sub>18</sub> F <sub>2</sub> N <sub>3</sub> O <sub>3</sub>
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

## Research

Area: Immunology SMILES: OC(=O)C1=CN(C2=CC=C(F)C=C2)C2=CC(N3CCNCC3)=C(F)C=C2C1=O Pathways: Others Receptor: Others Boiling pt: 621.4°C Melting pt: >240°C Solubility: DMSO: 1 mg/mL (2.37 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Cao Z, et al. Spectrochim Acta A Mol Biomol Spectrosc. 2015 Feb 5;136 Pt B:601-6.

Note: This product is for in vitro research use only