

2-Pyridinealdoxime methochloride

Catalog No: #S1111


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	2-Pyridinealdoxime methochloride
Brief Description	Inhibitors
Purification	98.00%
Target Name	AChR antagonist
Calculated MW	172.61
Formulation	C7H9ClN2O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
 SMILES: CN1C=CC=C/C1=C\[NH+]=O.[Cl-]
 Pathways: Neuroscience
 Receptor: AchE
 Boling pt: Melting pt: 230°C
 Solubility: DMSO: 10 mM
 Appearance: White Crystal
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Jokanović M, Prostran M. *Curr Med Chem.* 2009;16(17):2177-88.

Note: This product is for in vitro research use only