

## 2-(Hydroxymethyl)pyridine

Catalog No: #S1112


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	2-(Hydroxymethyl)pyridine
Brief Description	Inhibitors
Purification	98.00%
Target Name	PDE inhibitor
Calculated MW	109.13
Formulation	C6H7NO
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other  
SMILES: OCC1=CC=CC=N1  
Pathways: Metabolism  
Receptor: PDE  
Boiling pt: 112-113°C  
Melting pt: 5°C  
Solubility: DMSO: 10 mM  
Appearance: Liquid  
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

1. Ducharme Y, et al. Bioorg Med Chem Lett. 2003 Jun 2;13(11):1923-6.

Note: This product is for in vitro research use only