

Imazalil

Catalog No: #S1113


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Imazalil
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	297.19
Formulation	C ₁₄ H ₁₄ Cl ₂ N ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology
SMILES: ClC1=CC(Cl)=C(C=C1)C(CN1C=CN=C1)OCC=C
Pathways: Others
Receptor: Others
Boiling pt: >340°C
Melting pt: 52.7°C
Solubility: DMSO: > 10 mM
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Ruiz-Rodriguez L, et al. Food Chem. 2015 Jul 1;178:179-85.

Note: This product is for in vitro research use only