

Flumethasone

Catalog No: #S1124


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Flumethasone
Brief Description	Inhibitors
Purification	99.91%
Target Name	GR agonist
Calculated MW	410.46
Formulation	C22H28F2O5
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Endocrine

systemSMILES:C[C@@H]1C[C@H]2[C@@H]3C[C@@H](C4=CC(=O)C=C[C@@]4([C@]3([C@H](C[C@@]2([C@]1(C(=O)CO)O)C)O)F)C)FP
 Pathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:569.8C?Melting pt:237-240CSolubility:Ethanol: 6 mg/mL (14.61 mM); DMSO: 82 mg/mL (199.78 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Sternbauer K, et al. Zentralbl Veterinarmed A, 1998, 45(6-7), 441-443.

Note: This product is for in vitro research use only