

Penicillin V potassium salt

Catalog No: #S1126


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Penicillin V potassium salt
Brief Description	Inhibitors
Purification	99.95%
Target Name	Others
Calculated MW	388.48
Formulation	C16H17KN2O5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC1([C@@H](N2[C@H](S1)[C@@H](C2=O)NC(=O)COC3=CC=CC=C3)C(=O)[O-])C.[K+] Pathways: Others Receptor: Other
 sBoiling pt: 681.4°C Melting pt: 197-202°C Solubility: soluble in Water. Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Snchez-Morillas L, et al. Recent Pat Inflamm Allergy Drug Discov. 2014 Jan;8(1):74-6.

Note: This product is for in vitro research use only