

## Beclomethasone dipropionate

Catalog No: #S1127


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Beclomethasone dipropionate
Brief Description	Inhibitors
Purification	99.00%
Target Name	GR
Calculated MW	521.04
Formulation	C28H37ClO7
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Respiratory

systemSMILES:CCC(=O)OCC(=O)[C@]1([C@H](C[C@@H]2[C@@]1(C[C@@H]([C@]3([C@H]2CCCC4=CC(=O)C=C[C@@]43C)Cl)O)C)C)OC(=O)C

CPATHWAYS:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:630.5C?Melting pt:210CSolubility:DMSO: 10 mMAppearance:White

PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Saklatvala J. Arthritis Res. 2002;4(3):146-50.

Note: This product is for in vitro research use only