

Oxiracetam

Catalog No: #S1128


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Oxiracetam
Brief Description	Inhibitors
Purification	98.00%
Target Name	GABA Receptor agonist
Calculated MW	158.16
Formulation	C6H10N2O3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: NC(=O)CN1CC(O)CC1=O
Pathways: Membrane transporter/Ion channel;

Neuroscience
Receptor: GABAR
Boiling pt: 494.6°C
Melting pt: 165-168°C
Solubility: Water: 31 mg/mL (196.0 mM); DMSO: 31 mg/mL (196.0

mM)
Appearance: White Crystal
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.
Stock solution can be stored below -20 °C for several months.

References

1. Sansone M, et al. Arch Int Pharmacodyn Ther. 1985, 275(1), 86-92.

Note: This product is for in vitro research use only