

Nicergoline

Catalog No: #S1132


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nicergoline
Brief Description	Inhibitors
Purification	98.00%
Target Name	Adrenergic Receptor inhibitor
Calculated MW	484.39
Formulation	C ₂₄ H ₂₆ BrN ₃ O ₃
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Cardiovascular

System SMILES: CN1C[C@@H](C[C@]2([C@H]1CC3=CN(C4=CC=CC2=C34)C)OC)COC(=O)C5=CC(=CN=C5)Br

Pathways: Neuroscience; GPCR/G Protein Receptor: Adrenergic Receptor
Boiling pt: 594.4°C
Melting pt: Solubility: Soluble in DMSO
Appearance: White Powder
Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Alvarez-Guerra M, et al. Fundam Clin Pharmacol. 1999;13(1):50-8.

Note: This product is for in vitro research use only