

Clofoctol

Catalog No: #S1134


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Clofoctol
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	365.35
Formulation	C ₂₁ H ₂₆ Cl ₂ O
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Immunology SMILES: CC(C)(C)CC(C)(C)C1=CC(CC2=C(Cl)C=C(Cl)C=C2)=C(O)C=C1 Pathways: Others Receptor: Others Boiling pt: 447.7°C Melting pt: Solubility: DMSO: 10 mM Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

References

1. Yablonsky F, Simonnet G. J Pharmacol. 1982 Oct-Dec;13(4):515-24.

Note: This product is for in vitro research use only