

## Diastase

Catalog No: #S1136

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Diastase
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	N/A
Formulation	N/A
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Product Description

Research Area:Metabolism systemSMILES:N/APathways:OthersReceptor:StarchBoiling pt:Melting pt:Solubility:DMSO: 10 mMAppearance:Gray PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Gray GM. N Engl J Med. 1975 Jun 5;292(23):1225-30.

Note: This product is for in vitro research use only