

## Salbutamol


Catalog No: #S1139

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Salbutamol
Brief Description	Inhibitors
Purification	91.59%
Target Name	Adrenergic Receptor agonist
Calculated MW	239.31
Formulation	C13H21NO3
Storage	3 years -20°C powder;2 years -80°C in solvent;

## Images



## Product Description

Research Area:Respiratory systemSMILES:CC(C)(C)NCC(O)C1=CC(CO)=C(O)C=C1Pathways:GPCR/G ProteinReceptor:β-adrenergic receptorBoiling pt:433.5CMelting pt:157-158CSolubility:DMSO: 10 mMAppearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Bricchetto L, et al. Am J Physiol Lung Cell Mol Physiol. 2003 Jan;284(1):L133-9.

Note: This product is for in vitro research use only