

## Fenbendazole


Catalog No: #S1141

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Fenbendazole
Brief Description	Inhibitors
Purification	98.00%
Target Name	Others
Calculated MW	299.35
Formulation	C15H13N3O2S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: COC(=O)NC1=NC2=C(N1)C=C(SC1=CC=CC=C1)C=C2 Pathways: Others Receptor: Antiparasitic Boiling pt: Melting pt: 233°C Solubility: DMSO: 4 mg/mL (13.36 mM) Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Yu CG, et al. Neuroscience. 2014 Jan 3;256:163-9.

Note: This product is for in vitro research use only