

Dnmt3b Polyclonal Antibody

Catalog No: #40848

Package Size: #40848-1 50ul #40848-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Dnmt3b Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Hu
Specificity	Dnmt3b Polyclonal Antibody detects endogenous levels of Dnmt3b protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Dnmt3b.
Target Name	Dnmt3b
Other Names	DNMT3B; DNA (cytosine-5)-methyltransferase 3B; Dnmt3b; DNA methyltransferase HsallIB; DNA MTase HsallIB; M.HsallIB
Accession No.	Swiss-Prot: Q9UBC3NCBI Gene ID: 1789
Uniprot	Q9UBC3
GeneID	1789;
SDS-PAGE MW	96kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HT29 cells using Dnmt3b Polyclonal Antibody

Note: This product is for in vitro research use only