

## Diminazene Aceturate

Catalog No: #S1145


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Diminazene Aceturate
Brief Description	Inhibitors
Purification	95.00%
Target Name	Others; Serine Protease
Calculated MW	515.52
Formulation	C22H29N9O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Other SMILES: NC(=N)C1=CC=C(N=N)C2=CC=C(C=C2)C(N)=N)C=C1 Pathways: Others;

Proteases/Proteasome Receptor: Peroxiredoxin-5; QacR; Trypsin-1 Boiling pt: 463.3°C Melting pt: Solubility: Water: 100 mg/mL (193.97 mM); DMSO: 19 mg/mL (36.85 mM) Appearance: Yellow Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several months.

## References

- Overington JP, et al. Nat Rev Drug Discov. 2006 Dec;5(12):993-6.

Note: This product is for in vitro research use only