

Fusidine

Catalog No: #S1160


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Fusidine
Brief Description	Inhibitors
Purification	99.50%
Target Name	CAT; Others
Calculated MW	516.72
Formulation	C31H48O6
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research

Area: Immunology SMILES: CC(=O)OC12C[C@@H]3[C@H]([C@@H]1O)C[C@@H]2[C@@H]4[C@@H]3[C@@H]4C=C(C)C(=O)OC(=O)C Pathways: Others; oxidation-reduction Receptor: CAT; EEFGBoling pt: 635CMelting pt: 190-192CSolubility: Soluble in DMSO Appearance: White Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

References

1. Day PJ, et al. Biochemistry. 1995 May 16;34(19):6416-22.

Note: This product is for in vitro research use only