

Protionamide

Catalog No: #S1161


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-------------------|--|
| Product Name | Protionamide |
| Brief Description | Inhibitors |
| Purification | 97.30% |
| Target Name | Others |
| Calculated MW | 180.27 |
| Formulation | C9H12N2S |
| Storage | 3 years -20°C powder;2 years -80°C in solvent; |

Images



Product Description

Research Area:ImmunologySMILES:CCCC1=NC=CC(=C1)C(N)=SPathways:OthersReceptor:OthersBoiling pt:310.4CMelting pt:140-143CSolubility:Ethanol: 7 mg/mL (38.83 mM); DMSO: 36 mg/mL (199.7 mM)Appearance:Yellow PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Projahn M, et al. Antimicrob Agents ChemOthers. 2011 Sep;55(9):4408-11.

Note: This product is for in vitro research use only