

## Betamethasone 17-valerate

Catalog No: #S1163


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Betamethasone 17-valerate
Brief Description	Inhibitors
Purification	99.50%
Target Name	GR
Calculated MW	476.59
Formulation	C27H37FO6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:CCCCC(=O)O[C@@]1([C@H](C[C@@H]2[C@@]1(C[C@@H]([C@]3([C@H]2CCCC=CC(=O)C=C[C@@]43C)F)O)C)C(=O)COPat  
 hways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:598.9CMelting pt:183-184CSolubility:Ethanol: 58 mg/mL (121.7 mM); DMSO: 95 mg/mL (199.33 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Ortonne JP, et al. J Eur Acad Dermatol Venereol. 2014 Sep;28(9):1226-34.

Note: This product is for in vitro research use only