

## Desonide

Catalog No: #S1168


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Desonide
Brief Description	Inhibitors
Purification	99.50%
Target Name	GR antagonist
Calculated MW	416.51
Formulation	C24H32O6
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area:Endocrine

systemSMILES:C[C@]12C[C@@H]([C@H]3[C@H]([C@@H]1C[C@@H]4[C@]2(OC(O4)(C)C(=O)CO)CCC5=CC(=O)C=C[C@]35C)O)Pathways:MetabolismReceptor:Glucocorticoid ReceptorBoiling pt:Melting pt:272-274°C?Solubility:Ethanol: 30 mg/mL warmed (72.02 mM); DMSO: 83 mg/mL (199.27 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

## References

1. Russell GM, et al. J Neurosci. 2010 Apr 28;30(17):6106-15.

Note: This product is for in vitro research use only