

Moguisteine

Catalog No: #S1169


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Moguisteine
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others
Calculated MW	339.41
Formulation	C16H21NO5S
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:Respiratory systemSMILES:CCOC(=O)CC(=O)N1CCSC1COC1=CC=CC=C1OCPathways:OthersReceptor:OthersBoiling pt:506.8CMelting pt:72-76CSolubility:Ethanol: 68 mg/mL (200.34 mM); DMSO: 68 mg/mL (200.34 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

- Gallico L, et al. Br J Pharmacol, 1994, 112(3), 795-800.

Note: This product is for in vitro research use only