

## Diclazuril

Catalog No: #S1172


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Diclazuril
Brief Description	Inhibitors
Purification	98.77%
Target Name	Others
Calculated MW	407.64
Formulation	C17H9Cl3N4O2
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology  
 SMILES: ClC1=CC=C(C=C1)C(C#N)C1=C(Cl)C=C(C=C1Cl)N1N=CC(=O)NC1=O  
 Pathways: Others  
 Receptor: Others  
 Boiling pt: Melting pt: 548°C  
 Solubility: DMSO: 16 mg/mL (39.25 mM)  
 Appearance: White Powder  
 Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Oz HS, et al. Int J Clin Med. 2014 Jan 1;5(3):93-101.

Note: This product is for in vitro research use only