

Racecadotril


Catalog No: #S1176

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Racecadotril
Brief Description	Inhibitors
Purification	99.39%
Target Name	Proteasome inhibitor
Calculated MW	385.49
Formulation	C ₂₁ H ₂₃ NO ₄ S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Metabolism system
SMILES: CC(=O)SCC(CC1=CC=CC=C1)C(=O)NCC(=O)OCC1=CC=CC=C1
Pathways: Proteases/Proteasome; Ubiquitination
Receptor: Endopeptidase
Boiling pt: 574.5°C
Melting pt: 89°C
Solubility: DMSO: 22mg/ml (57.07 mM); Water:

References

1. Matheson AJ, Noble S. *Drugs*. 2000 Apr;59(4):829-35; discussion 836-7.

Note: This product is for in vitro research use only