

## Sulfaquinoxaline

Catalog No: #S1177


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Sulfaquinoxaline
Brief Description	Inhibitors
Purification	98.35%
Target Name	Others
Calculated MW	300.34
Formulation	C <sub>14</sub> H <sub>12</sub> N <sub>4</sub> O <sub>2</sub> S
Storage	3 years -20°C powder; 2 years -80°C in solvent;

## Images



## Product Description

Research Area: Immunology SMILES: NC1=CC=C(G=C1)S(=O)(=O)NC1=NC2=CC=CC=C2N=C1 Pathways: Others Receptor: Others Boiling pt: Melting pt: 248-255°C Solubility: Water: 7.5 mg/L Appearance: Khaki Powder Remark: For obtaining a higher solubility, please warm the tube at 37 °C and shake it in the ultrasonic bath for a while. Stock solution can be stored below -20 °C for several month.

## References

1. Guo Y, et al.  
Anal Chem. 2010 Sep 15;82(18):7550-5.

Note: This product is for in vitro research use only