

Blonanserin

Catalog No: #S1180


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Blonanserin
Brief Description	Inhibitors
Purification	99.50%
Target Name	5-HT Receptor antagonist; Dopamine Receptor antagonist
Calculated MW	367.5
Formulation	C23H30FN3
Storage	3 years -20°C powder; 2 years -80°C in solvent;

Images



Product Description

Research Area: Nervous system
SMILES: CCN1CCN(CC1)C1=NC2=C(CCCCC2)C(=C1)C1=CC=C(F)C=C1
Pathways: GPCR/G Protein;
Neuroscience Receptor: 5-HT₂; D₂
Boiling pt: 540.8°C
Melting pt: 117-119°C
Solubility: Ethanol:

References

1. Kawabe K, et al. Clin Neuropharmacol. 2013, Nov-Dec; 36(6):239-41.

Note: This product is for in vitro research use only