

Bifonazole

Catalog No: #S1187


Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Bifonazole
Brief Description	Inhibitors
Purification	99.50%
Target Name	Others; P450 inhibitor
Calculated MW	310.39
Formulation	C22H18N2
Storage	3 years -20°C powder;2 years -80°C in solvent;

Images



Product Description

Research Area:ImmunologySMILES:C1=CN(C=N1)C(C1=CC=CC=C1)C1=CC=C(C=C1)C1=CC=CC=C1Pathways:Metabolism;
OthersReceptor:Others; CYP51Boiling pt:491.7CMelting pt:142CSolubility:Ethanol: 20 mg/mL (64.43 mM); DMSO: 62 mg/mL (199.74 mM)Appearance:White PowderRemark:For obtaining a higher solubility , please warm the tube at 37 °C and shake it in the ultrasonic bath for a while.Stock solution can be stored below -20 °C for several month.

References

1. Carrillo-Munoz AJ, et al. Rev Esp Quimioter. 2006 Jun;19(2):130-9.

Note: This product is for in vitro research use only